



Prosthetic Lamination Lay-ups

Lay Up 1

**½ oz Dacron
1 Fiberglass Tubing
4 Nyglass
1 Fiberglass Tubing
2 Nyglass
Finish Nylon**

Lay Up 2

**½ oz Dacron
1 Fiberglass Tubing
6 Nyglass
1 Fiberglass Tubing
2 Nyglass
Finish Nylon**

Lay Up 3

**½ oz Dacron
1 Fiberglass Tubing
10 Nyglass
1 Fiberglass Tubing
2 Nyglass
Finish Nylon**

Lay Up 4

**½ oz Dacron
1 Carbon Tubing
4 Nyglass
1 Carbon Tubing
Finish Nylon**

Lay Up 5

**CarbonHvy
½ oz Dacron
1 Carbon Tubing
6 Nyglass
1 Carbon Tubing
Finish Nylon**

Lay Up 6

**AK Frame Lay UP
Lam 1
½ oz Dacron
1 Fiberglass Tubing
4 Nyglass**

Lam 2

**1 Fiberglass Tubing
4 Nyglass
Finish Nylon**